



REALISTIC. DYNAMIC.
Organ-on-Chip Models



SynBBB


Blood-Brain-Barrier Model

TEER CONFIGURATION

For Chips and Starter Kits:
Catalog #s: 402004, 102015-SB3

INTRODUCTION	3
SynVivo Cell Impedance Analyzer	4
Cell Impedance Analyzer Specifications	4
Choosing An Operating Frequency	4
SynVivo SynBBB Protocols	5-17
Workflow	5
Basement Membrane	6-8
Priming SynBBB	9-10
Z Measurement	11
Establishing Monoculture	12-13
Establishing Tri-Culture	14-15
3D Shear Monoculture Model	16
3D Shear Tri-Culture Model	17
Additional Resources – Peer-reviewed Publications	18

Introduction



Schematic of the chip used for the SynBBB Model and SynVivo Cell Impedance Analyzer system. Vascular chambers (outer channels) are for culture of endothelial cells while the tissue chamber (central chamber) is for culture of brain tissue cells (astrocytes, pericytes, neurons). Porous architecture enables communication between the vascular and tissue cells.

Delivery of neuroprotective or therapeutic agents to specific regions of the brain presents a major challenge, largely due to the presence of the Blood-Brain Barrier (BBB). Physiologically, the BBB consists of an intricate network of vascular endothelial cells that isolate the central nervous system from systemic blood circulation. Traditional BBB assays, such as the 2D cell culture insert model have significant limitations such as lack of physiological shear stress, real-time visualization capability and large number of consumables. SynVivo SynBBB recreates the in vivo microenvironment by mimicking a histological slice of brain tissue cells in communication with endothelial cells across the BBB. SynBBB is the only in vitro BBB model that allows:

- Accurate in vivo hemodynamic shear stress
- Real-time visualization of cellular and barrier functionality
- Significant reduction in cost and time
- Robust and easy to use protocols

This model has been successfully demonstrated for upregulation of tight junction molecules, functional assays (Prabhakarapandian et al., 2013) and validated against in vivo studies showing excellent correlation with permeation of small molecules (Deosarkar et al., 2015). Electrical resistance measurements provide a non-invasive method for real-time monitoring of tight junctions as detailed below. Tight junction formation between cells (e.g., blood brain barrier) can be evaluated by measuring changes in electrical resistance in the intercellular space between cells. The SynVivo Cell Impedance Analyzer, used in conjunction with the SynBBB TEER chip, measures electrical impedance (resistance, denoted by "Z").

SynVivo Cell Impedance Analyzer Specifications

Impedance Measurement Range


Minimum	20 kOhm
Maximum	300 kOhm

Operating Frequencies

500, 1,000 and 10,000 Hz

Connections

Electrode	Micro-grabber Clips
Power	AC Adapter



The SynVivo Cell Impedance Analyzer is designed to measure changes in impedance-base resistance (Z) across cells cultured in SynVivo SynBBB chips. A toggle switch on the right side of the Analyzer changes the operational frequency between 500, 1,000, and 10,000 Hz. The Analyzer is interfaced to the SynVivo SynBBB chip by attaching micro-grabber clips to AgCl electrodes inserted into the electrode ports of the chip. The measurement button is pressed once to display the operational frequency, and once again to acquire the resistance Z and phase measurement.

Important Note: The SynVivo® Cell Impedance Analyzer is compatible only with SynBBB TEER chips.

Choosing an Operating Frequency

The table below gives optimized operating frequencies for certain cell types.

Cell Type	Operating Frequency
Human Umbilical Vein Endothelial Cell (HUVEC)	10 kHz
Human Brain Microvascular Endothelial Cells	1 kHz
Rat Brain Microvascular Endothelial Cells (RBEC)	1 kHz


Use the following protocol to determine the optimal operating frequency for a specific cell type:

1. Measure the baseline impedance and phase values of a coated and fully perfused chip at operating frequencies of 500, 1,000, and 10,000 Hz.
2. Coat the chip with the desired cells and grow under desired conditions.
3. Acquire impedance and phase values at operating frequencies of 500, 1,000, and 10,000 Hz.
4. Choose the operating frequency with a phase value closest to 0 degrees, ideally between 0 and ± 15 degrees.


SynVivo SynBBB Protocols

Workflow

The following workflow charts shows the timing of establishing monoculture and tri-culture SynBBB TEER models:



Monoculture Workflow



Tri-culture Workflow

Establishing the Basement Membrane

Principle:


The extracellular matrix (ECM), or basement membrane, is a structural framework of proteins that provides a scaffold for cells comprising a tissue. Fibronectin is a major component of vascular endothelial networks, where it serves as a scaffold protein and plays a role in vascular morphogenesis. This protocol will use human fibronectin as the ECM, but the protocol can be adapted for other types of ECM, such as collagen or Matrigel.

Equipment:

- SynVivo SynBBB Impedance (TEER) Chips (SynVivo 102015-SB3)
- Tygon tubing (0.02 inch ID x 0.06 inch OD, or 0.05 cm ID x 0.15 cm OD; SynVivo 201005)
- Blue slide clamps (SynVivo 202001)
- 1 mL syringe with Luer-Lock tip (25 pack; SynVivo 203005)
- 24-gauge blunt-tipped needles (0.5 inches or 1.27 cm long; SynVivo 204003)
- Pneumatic Primer (SynVivo 205001)
- Forceps
- Scissors


Reagents:

- 1X PBS without calcium or magnesium
- Serum-free endothelial cell basal medium
- Human fibronectin



Protocol for Monoculture:

1. Using forceps, insert 1-inch segments of Tygon tubing into the outlet **Ports (B, D, and F)** and the electrode **Ports (1, 2, and 3)**. Serum-free endothelial cell basal medium
2. Thaw and dilute fibronectin to 200 $\mu\text{g}/\text{mL}$ in chilled serum-free media or 1X PBS.
3. Attach a 24-gauge blunt-tipped needle to a 1ml syringe. Cut and attach an approximately 5-inch segment of Tygon tubing to the end of the needle.
4. Draw up the diluted fibronectin only into the Tygon tubing (not into the needle or syringe). Attach the free end of the tubing to inlet **Port A** of the vascular channel.
5. Gently push fibronectin into the vascular channel until 2 droplets form on the outlet tubing in **Port B**. Cut the tubing connected to the syringe approximately 1 inch from the chip.
6. Repeat this process for the other vascular channel, using inlet **Port E** and outlet **Port F**.
7. Place the chip at 37°C in an incubator for 1 hour. The chip can go as long as overnight, if needed.
8. The chip is ready to be primed with the Pneumatic Primer.



Protocol for Tri-culture:

1. Using forceps, insert 1-inch segments of Tygon tubing into the outlet **Ports (B, D, and F)** for the vascular channel and the electrode **Ports (1, 2, and 3)**. Clamp the tubing attached to electrode **Port 2**.
2. Thaw and dilute fibronectin to 200 $\mu\text{g}/\text{mL}$ in chilled serum-free media or 1X PBS.
3. Attach a 24-gauge blunt-tipped needle to a 1ml syringe. Cut and attach an approximately 5-inch segment of Tygon tubing to the end of the needle.
4. Draw up the diluted fibronectin only into the Tygon tubing (not into the needle or syringe). Attach the free end of the tubing to inlet **Port A** of the vascular channel.
5. Gently push fibronectin into the vascular channel until 2 droplets form on the outlet tubing in **Port B**. Cut the tubing connected to the syringe approximately 1 inch from the chip.
6. Repeat this process for the other vascular channel (inlet **Port E** and outlet **Port F**) and the tissue channel (inlet **Port C** and outlet **Port D**).
7. Place the chip at 37°C in an incubator for 1 hour. The chip can go as long as overnight, if needed.
8. The chip is ready to be primed with the Pneumatic Primer.

Priming SynBBB TEER Chips

Principle:




Complete removal of air from the SynBBB chips is critical for both the successful culture of cells and accurate measurement of resistance using the Cell Impedance Analyzer. The Pneumatic Primer (SynVivo 205001) works by pressurizing the SynBBB chip and forcing trapped air bubbles through the pores in the PDMS. This method requires a source of inert gas, as well as a pressure regulator to control the flow of gas into the chip.

Equipment:

- SynVivo SynBBB Impedance (TEER) Chips (SynVivo 102015-SB), coated with ECM of choice
- Tygon tubing (0.02 inch ID x 0.06 inch OD, or 0.05 cm ID x 0.15 cm OD; SynVivo 201005)
- Blue slide clamps (SynVivo 202001)
- 1 mL syringe with Luer-Lock tip (25 pack; SynVivo 203005)
- 24-gauge blunt-tipped needles (0.5 inches or 1.27 cm long; SynVivo 204003)
- Pneumatic Primer (SynVivo 205001)
- Forceps
- Scissors
- Multiple port manifold for the Pneumatic Primer (SynVivo 20700)
- Source of inert compressed gas

Reagents:

- Serum-free endothelial cell basal medium



Protocol:

1. Attach a 24-gauge blunt-tipped needle to a 1ml syringe. Fill the syringe with either serum-free basal medium.
2. Cut and attach a 5-8 inch segment of Tygon tubing to the end of the needle, and using gentle pressure, fully perfuse the attached tubing.
3. Connect the tubing to **Port A**. Gently push liquid into the chip until 1-2 drops of liquid come out the tubing of **Port B**. Clamp the **Port B** tubing below the liquid line. Continue to gently push liquid into the chip until tubing attached to electrode **Port 1** is filled with liquid. Clamp **Port 1** tubing below the liquid line.
4. Cut the tubing attached to the syringe to approximately 1-inch from **Port A** and clamp.
5. Connect the tubing to **Port E**. Gently push liquid into the chip until 1-2 drops of liquid come out the tubing of **Port F**. Clamp the **Port B** tubing below the liquid line. Continue to gently push liquid into the chip until tubing attached to electrode **Port 3** is filled with liquid. Clamp **Port 3** tubing below the liquid line.
6. Cut the tubing attached to the syringe to approximately 1-inch from **Port E** and clamp.
7. Connect the tubing to **Port C**. Gently push liquid into the chip until 1-2 drops of liquid come out the tubing of electrode **Port 2**. Clamp the **Port 2** tubing below the liquid line. Continue to gently push liquid into the chip until tubing attached to **Port D** is filled with liquid. Clamp **Port D** tubing below the liquid line.
8. With one end of Tygon tubing still attached to **Port C**, slip the other end of the tubing off the needle. All port tubes except **Port C** should be clamped.
9. Connect the chip to the Pneumatic Primer by attaching the **Port C** tubing to the needle on the Luer-Lock connector on the manifold. Note: Multiple chips can be primed simultaneously using the multiple port manifold, available from SynVivo (cat # 207001)
10. Turn the knob on the controller box and adjust the pressure to 7.5 psi. Apply the pressure for ~20-25 minutes.
11. Turn off the pressure and cut the Tygon tubing connected to the Pneumatic Primer.
12. Allow the chip to incubate at 37°C for a minimum of 20 minutes before introducing cells.

Measuring Resistance Z

Principle:

The SynVivo Cell Impedance Analyzer is used measure changes in electrical resistance between cell layers as a method of monitoring cellular tight junctions. To measure the resistance of the cell layer, acquire measurements before and after culturing cells in the primed SynBBB chip. Calculated the cell layer Z by subtracting the baseline resistance of the primed SynBBB chip (no cells) from the resistance of a SynBBB chip functionalized with cells. This is equivalent to subtracting the TEER of a blank cell culture insert from TEER of a functionalized insert in a 2-D cell culture system.

Equipment:

- SynBBB Impedance (TEER) Chips (Catalog # 102015-SB)
- SynVivo Impedance Analyzer (Catalog # 304001)
- AgCl Electrodes
- Hotplate (optional)

Protocol:

1. Place a fully perfused SynBBB chip on a hotplate set to 37°C. If a hotplate is not available, allow the chip to come to room temperature for ~10 minutes before proceeding.
2. Clean two AgCl electrodes with 70% ethanol or isopropanol.
3. Using forceps, thread an electrode through electrode **Port 2** on the chip, pushing the electrode tip down in contact with the glass slide.
4. Thread the 2nd electrode into electrode **Port 1**. Ensure that both electrode tips are touching the glass slide.
5. Turn on the SynVivo Cell Impedance Analyzer by plugging in the power cord. The Analyzer will calibrate automatically.
6. Attach a micro-grabber clip to each electrode, being sure not to pull the electrodes out of place.
7. Select the desired operating frequency using the toggle located on the right side of the Analyzer (**see Choosing an Operating Frequency on Page 4**). For primary human brain microvascular endothelial cells, choose 1,000 Hz. Press the measurement button to display the impedance (Z) and phase values on the screen.
8. If another frequency is desired, switch the toggle and press the measurement button, once to display the current operating frequency, and once again to display the impedance and phase values.
9. Repeat steps 3 through 7 on electrode **Port 3**.

Establishing Monoculture

Principle:

Shear-induced endothelial cell tight junctions, which cannot be achieved in a 2-D model, are readily achieved in the SynBBB model using fluid perfusion. Formation of tight junction changes can be measured using biochemical or electrical analysis (assessing changes in electrical resistance Z) with the SynVivo Cell Impedance Analyzer (SynVivo 304001).

Note: *Primary human brain microvascular endothelial cells (HBMVECs) are subject to phenotypic drift with subsequent passaging. To preserve the phenotypic tight junctions, do not passage beyond the manufacturer's instructions.*

Equipment:

- SynVivo Pneumatic Primer Device (Catalog #205001)
- 1 mL BD plastic syringes or other 1 mL syringes (Catalog # 203004)
- 24-gauge blunt tip needles (Catalog # 204002)
- Tygon microbore tubing, 0.02" ID X 0.06" OD (Catalog # 201005)
- Clamps (Catalog # 202003)
- Forceps
- Syringe Pump capable of flow rates from 10nl/min to 10μlmin
- Hemacytometer

Reagents:

- Primary human brain microvascular endothelial cells
- Complete Endothelial Media
- Dissociation reagents and neutralizers
- 1X PBS without calcium or magnesium

Establishing Monoculture

Protocol:

1. Endothelial cells should be dissociated, centrifuged, and concentrated to approximately 3×10^7 cells/mL. **Note:** Do not expose primary cells to trypsin longer than necessary. Most primary cells will dissociate in <3 minutes.
2. Place a drop of PBS at the base of the inlet **Port A** tubing to be removed.
3. Gently remove the tubing—the liquid droplet should cover the port once the tubing is removed, and should remain there until new Tygon tubing is inserted into the port. This will prevent air from entering the chip.
4. Remove the clamp on the outlet **Port B**.
5. Attach a 24-gauge blunt-tip needle to a 1 mL syringe. Cut 5-8 inches of Tygon tubing and attach one end to the needle.
6. Draw the cells up only into the Tygon tubing and mount the syringe onto a syringe pump.
7. Ensure that the tubing is free of air bubbles and allow a convex meniscus to form at the end of the tubing.
8. Insert the tubing into the **Port A**—the meniscus and the drop of water will form a bridge to prevent air from entering the chip as the tubing is inserted.
9. Clean the fluid from the surface of the chip.
10. Begin the injection at 3 $\mu\text{L}/\text{min}$.
11. Observe the chip as the cells are flowing. Once the vascular channel is filled with cells, clamp the outlet tubing (**Port B**) and stop the flow.
12. Cut the inlet tubing to separate the chip from the syringe.
13. Repeat this process for the other vascular channels using **Port E and F**. Cells should be uniformly dispersed at this moment.
14. Place the chip in the incubator for a minimum of 4 hours before inducing continuous flow.
15. Continuous flow of cell media is slowly ramped via programmable syringe pump overnight using the following steps:
 - a) Prepare a 1 mL syringe with endothelial cell media and attach a 24-gauge blunt-tip needle and a 12 to 20-inch length of Tygon tubing.
 - b) Mount the syringe to a programmable syringe pump and connect the Tygon tubing to the SynBBB chip in the incubator.
 - c) Program the following:
 - i. Step 1: infuse 0.02 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - ii. Step 2: infuse 0.03 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - iii. Step 3: infuse 0.04 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - iv. Step 4: infuse 0.05 $\mu\text{L}/\text{min}$, 14:00:00 (14 hours)
 - d) Start the flow program.
16. For a tri-culture model, continue with “Establishing Tri-Culture”. For a monoculture model, skip to “3D Shear Monoculture Model”.

Establishing Tri-culture

Principle:

HBMVECs cultured alongside astrocytes and pericytes will form tighter junctions than brain endothelial cells alone. The following protocol is an example of tri-culture using primary brain endothelial cells cultured alongside primary human astrocytes and pericytes. This protocol may be adapted for any vascular-tissue model.

Equipment:

- SynVivo Pneumatic Primer Device (Catalog #205001)
- 1 mL BD plastic syringes or other 1 mL syringes (Catalog # 203004)
- 24-gauge blunt tip needles (Catalog # 204002)
- Tygon microbore tubing, 0.02" ID X 0.06" OD (Catalog # 201005)
- Clamps (Catalog # 202003)
- Forceps
- Syringe Pump capable of flow rates from 10nl/min to 10µl/min
- Hemacytometer

Reagents:

- Human Astrocytes
- Human Pericytes
- Complete Astrocyte Media
- Complete Pericyte Media
- Dissociation reagents and neutralizers
- 1X PBS without calcium or magnesium or
- Serum-free basal medium

Establishing Tri-culture

Protocol:

1. All tubing should be clamped before beginning this procedure.
2. Dissociate and re-suspend astrocytes and pericytes at a concentration of approximately 1×10^7 cells/mL. Mix the two cell types in a 2:1 ratio of astrocytes:pericytes.
3. Place a drop of 1X PBS at the base of inlet **Port C**.
4. Gently remove the tubing - the PBS should cover the port once the tubing is removed and should remain there until new Tygon tubing is inserted into the port to prevent air from entering the chip.
5. Remove the clamp on the outlet **Port D**.
6. Prepare a syringe and tubing with the cell suspension and mount onto a syringe pump.
7. Push the cell suspension to the end of the tubing to form convex meniscus.
8. Insert the tubing into the **Port C**—the meniscus and the drop of water will form a bridge to prevent air from entering the chip as the tubing is inserted.
9. Begin the injection at 3 μ L/min.
10. Observe the chip as the cells are flowing. Once the tissue chamber is filled with cells, stop the flow and clamp the outlet tubing (**Port D**).
11. Cut the inlet tubing and allow the chip to incubate at 37°C for **4 hours** before resuming flow (**3D Shear Tri-culture Model**).

3D Shear Monoculture Model

Principle:

Flow is continued on monoculture models 24 hours after seeding the primary brain endothelial cells. The monoculture model is ready for use 48 hours after the initial seeding of endothelial cells.

Equipment:

- SynBBB chip with endothelial cells
 - Programmable syringe pump
 - 1 mL syringe
 - 24-Gauge blunt-tip needle
 - Tygon tubing Reagents:
 - Endothelial Cell Culture Media
 - 1X PBS
1. Prepare a 1 mL syringe with endothelial cell media and attach a 24-gauge blunt-tip needle and a 12 to 20-inch length of Tygon tubing.
 2. Mount the syringe to a programmable syringe pump and connect the Tygon tubing to the SynBBB chip in the incubator.
 3. Program the following:
 - a) **Step 1:** infuse 0.05 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - b) **Step 2:** infuse 0.06 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - c) **Step 3:** infuse 0.07 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - d) **Step 4:** infuse 0.08 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - e) **Step 5:** infuse 0.08 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - f) **Step 6:** infuse 0.1 $\mu\text{L}/\text{min}$, 14:00:00 (14 hours)
 4. The SynBBB model is ready for assay and resistance Z measurement (**see “Measuring Resistance Z” on page 11**).

3D Shear Tri-Culture Model

Principle:

Flow is continued on tri-culture models 4 hours after seeding astrocytes and pericytes in the tissue chamber. The tri-culture model is ready for use 72 hours after the initial seeding of endothelial cells.

Equipment:

- SynBBB chip with endothelial cells
 - Programmable syringe pump
 - 1 mL syringe
 - 24-Gauge blunt-tip needle
 - Tygon tubing Reagents:
 - Endothelial Cell Culture Media
 - 1X PBS
1. Prepare a 1 mL syringe with endothelial cell media and attach a 24-gauge blunt-tip needle and a 12 to 20-inch length of Tygon tubing.
 2. Mount the syringe to a programmable syringe pump and connect the Tygon tubing to the SynBBB chip in the incubator.
 3. Program the following:
 - a) **Step 1:** infuse 0.05 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - b) **Step 2:** infuse 0.06 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - c) **Step 3:** infuse 0.07 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - d) **Step 4:** infuse 0.08 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - e) **Step 5:** infuse 0.08 $\mu\text{L}/\text{min}$, 2:00:00 (2 hours)
 - f) **Step 6:** infuse 0.1 $\mu\text{L}/\text{min}$, 34:00:00 (34 hours)
 4. The SynBBB model is ready for assay and resistance Z measurement (**see “Measuring Resistance Z” on page 11**).

Additional Resources – Peer-Reviewed Publications

For further information on how the SynBBB model has been used in research, see the following supporting information:

1. Greene C, Rebergue N, Fewell G, Damir Janigro, Yann Godfrin, Campbell M, Sighild Lemarchant. 2024. NX210c drug candidate peptide strengthens mouse and human blood-brain barriers. *Fluids and Barriers of the CNS*. 21(1). doi:<https://doi.org/10.1186/s12987-024-00577-x>.
2. Huang J, Li YB, Charlebois C, Nguyen T, Liu Z, Bloembergen D, Zafer A, Baumann E, Sodja C, Leclerc S, et al. 2022. Application of blood brain barrier models in pre-clinical assessment of glioblastoma-targeting CAR-T based immunotherapies. *Fluids and Barriers of the CNS*. 19(1). doi:<https://doi.org/10.1186/s12987-022-00342-y>.
3. Brown TS, Nowak M, Bayles AV, Balabhaskar Prabhakarapandian, Pankaj Karande, Lahann J, Helgeson ME, Samir Mitragotri. 2019. A microfluidic model of human brain (μ HuB) for assessment of blood brain barrier. *Bioengineering & translational medicine*. 4(2). doi:<https://doi.org/10.1002/btm2.10126>.
4. Tang Y, Soroush F, Sun S, Liverani E, Langston JC, Yang Q, Kilpatrick LE, Kiani MF. 2018. Protein kinase C-delta inhibition protects blood-brain barrier from sepsis-induced vascular damage. *Journal of Neuroinflammation*. 15(1). doi:<https://doi.org/10.1186/s12974-018-1342-y>.
5. Terrell-Hall TB, Nounou MI, Fatema El-Amrawy, Griffith JIG, Lockman PR. 2017. Trastuzumab distribution in an in-vivo and in-vitro model of brain metastases of breast cancer. *Oncotarget*. 8(48):83734–83744. doi:<https://doi.org/10.18632/oncotarget.19634>.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses.