







# Assay Kit Contents

This list applies to the following catalog numbers: 401001 \* 401003 \* 402001 \* 402003 \* 403001 \* 403003 \* 403005 \* 403007 \* 404001

<p>(A) SynVivo Chips (10)</p> <p>Cat #: variable</p>		<p>The chips are individually packaged in resealable silver bags.</p>
<p>(B) Tubing (100 ft)</p> <p>Cat #: 201005</p>		<p>Tubing will connect the ports of the microchip to pumps, pneumatic primers, etc.</p>
<p>(C) Slide Clamps (25)</p> <p>Cat #: 202003</p>		<p>25 blue, rectangular slide clamps will be included, and are intended to seal the microchip system, for example, after cells have been seeded.</p>
<p>(D) Blunt Tip Needles (50)</p> <p>Cat #: 204002</p>		<p>50 needles are included to connect the syringes to the tubing.</p>




(E) 1 mL  
Syringes  
(50)


Cat #:  
203004



50 syringes are included to enable various media flows through a “feeding” pump.



Prime



Seeding Closeup



Media Flow

MANIFOLD, PNEUMATIC PRIMER AND ADAPTER NOT INCLUDED